

PRODUCT DISRUPTION BRIEF

Tourniquet cuffs (updated December 13)

Tourniquet cuff market has recovered, all products are in stock and ready to ship.

Current conditions

Vizient has confirmed that Stryker is currently shipping and stocking all products previously on backorder or in safety stock:

Size	Part Numbers	Status
9"	5921-209-235, 5921-209-135	Ready to ship
15"	5921-115-135, 5921-115-235	Ready to ship
18"	5921-018-245, 5921-218-235, 5921-018-145, 5921-018-235	Ready to ship
24"	5921-024-245, 5921-024-235, 5921-024-145, 5921-024-135	Ready to ship
30"	5921-030-235, 5921-030-135	Ready to ship
34"	5921-034-235, 5921-034-135	Ready to ship

Coming Soon SKUs:

Size	Part Numbers	Status
Reusable	all sizes	Backorder cleared and ready to ship
9"	5921-109-135, 5921-109-235	Backorder cleared and ready to ship
12"	5921-212-235, 5921-212-135, 5921-112-135, 5921-112-235	In safety stock and ready to ship
44"	5922-044-135, 5922-044-235, 5921-018-145, 5921-018-245, 5921-024-145, 5921-024-245	Backorder cleared and ready to ship

Vizient has completed the bid process and the following suppliers have been awarded, effective Feb 2024: ASP Global (MS1272), SourceMark (MS1273), Stryker (MS1274) and Zimmer (MS1275).

In February 2023, Stryker issued an [urgent medical device recall](#) for specific lots of color cuff non-sterile disposable tourniquet cuffs. Reported incidences indicate reporting of low-pressure events including not holding pressure and not achieving desired pressure during initial set up. This is caused by the fabric of the cuff tearing during set up and a pinched or unseated O-ring at the connection between the cuff, disposable fill line and/or cuff housing. Stryker has contacted all affected members. Healthcare providers should review the Safety Directives and Instructions For Use (IFU) available in the [urgent medical device recall notice](#) to mitigate the risk to achieve or maintain pressure. If product does not reach or maintain targeted pressure; discontinue use and refer to action steps provided in the [urgent medical device recall notice](#). If product does maintain pressure, healthcare providers can continue to use or request a credit.

Alternatives

Refer cross references of Stryker at the following links: [non-sterile products](#) and [sterile products](#). A downloadable Excel version is available [here](#).

Alternative supply may be sourced from the following off-contract suppliers; Additional commitment may be required.

	Delfi
Product name	Personalized Tourniquet System (PTS)
Offering	tourniquet pump, cuffs (single, dual, reusable), sleeves (stockinette, sizes (3"-52"))
Website	Home - Delfi (delfimedical.com)
Customer svc	1-800-933-3022, request sales at greeting
Info	Delfi info sheet

Compatibility/adapters

Stryker has provided additional information regarding cuff compatibility with Zimmer tourniquet systems [here](#). All Stryker cuffs produced are tested and safe for use with the Zimmer tourniquet systems (Zimmer tourniquet pumps – Models ATS500,



To learn more, please contact:
Hadeer Fati, Category Mgr.
hadeer.fati@vizientinc.com.

As the nation's largest member-driven health care performance improvement company, Vizient provides solutions and services that empower health care providers to deliver high-value care by aligning cost, quality and market performance. With analytics, advisory services and a robust sourcing portfolio, we help members improve patient outcomes and lower costs.

ATS750, ATS1200, ATS2000, ATS3000, and ATS4000). These cuffs have the same proprieties and functions as any other tourniquet cuff on the market and will perform within all clinical performance specifications for the established cuffs.

Stryker also offers the following adapters with healthy available stock; work with your local sales representative to determine which adapter is most appropriate for your situation.

The below adapters are designed to connect to single port cuffs:

- 5920-000-301 Adapter, single, quick lock to luer lock
- 5920-000-305 Adapter set, quick lock to luer lock

The below adapters are designed to connect to dual port cuffs:

- 5920-000-302 Y adapter, single, luer lock to luer lock
- 5920-000-303 Y adapter, single, quick lock to luer lock
- 5920-000-304 Y adapter, single, quick lock to quick lock
- 5920-000-306 Y adapter set, quick lock to luer lock
- 5920-000-307 Y adapter set, quick lock to quick lock

Mitigation strategies

In June, the FDA released the following [Non-sterile, Single-use Pneumatic Tourniquet Cuffs Conservation Strategies - Letter to Health Care Providers](#) to address the use of alternative devices and existing devices after appropriate sterilization/disinfection procedures have been appropriately followed.

According to the articles linked in this paragraph, a single-use sterile [Class 1 medical device](#), registered with the FDA can be used for limb exsanguination during limb surgery to provide a bloodless surgical field. The device consists of an elastic silicone ring, a stockinette and straps with handles. By pulling on the straps, it can be unrolled along the limb proximal to the surgical area. It can be used for all surgical procedures involving limbs where a bloodless field is needed, including trauma, orthopaedic, hand, foot and ankle, plastic and vascular cases. [here](#), [here](#) and [here](#).

Reprocessing capability

Collections are key to any reprocessing program. Start collecting any single use tourniquet cuffs immediately instead of discarding. Product back to a Member should be expected only for that which has been collected and will only begin after product has been collected.

Instructions for collections must be communicated with Member prior to the go-live of any new reprocessing program to ensure success of program and shipment of product back to Member.

- **Stryker Sustainability** offers reprocessed products for all Stryker and Zimmer product codes including previous Stryker exceptions: 5921-030-135, 5921-030-235, 5921-034-135, 5921-034-235, 5921-218-135 and 5921-218-235. OEM Stryker cuffs can be reprocessed up to 3 times; OEM Zimmer cuffs are cleared to be reprocessed up to 6 times. Current customers should work with your current sustainability rep. Members not currently reprocessing should reach out to their Nathan.Javadi@stryker.com or their local Stryker rep for information on how to start collections & reprocessing.
- **Medline** has the capability to reprocess Stryker and Zimmer tourniquet cuffs. Medline's clearance has been approved to allow reprocessing up to 5 times for both Zimmer & Stryker cuffs. Members not currently reprocessing with Medline should contact Jake Colletti @ JColletti@medline.com.
- **Cardinal** can reprocess Stryker tourniquet cuffs; however, members should work with Cardinal to confirm SKU allowances for reprocessing. Cardinal has no apparent capabilities for reprocessing Zimmer tourniquet cuffs. Members should discuss with their local Cardinal rep for more information and to connect with a reprocessing specialist.
- **Arjo** has the ability to reprocess only some Stryker tourniquet cuffs. Members should connect with their local Arjo rep for additional information.

DISCLAIMER: THE INFORMATION CONTAINED IN THIS DOCUMENT IS INTENDED FOR INFORMATIONAL PURPOSES ONLY AND IS IN NO WAY INTENDED TO BE A SUBSTITUTE FOR OR IN ANY MANNER TO BE CONSTRUED AS MEDICAL OR CLINICAL ADVICE. VIZIENT IS COMPILING INFORMATION AND EMERGING PRACTICES FROM MEMBERS TO AID IN KNOWLEDGE TRANSFER DURING CRITICAL SUPPLY EVENTS. THE INFORMATION CONTAINED HEREIN HAS NOT BEEN INDEPENDENTLY VERIFIED, RESEARCHED, OR INVESTIGATED AND SHOULD NOT BE CONSTRUED AS ADVICE OR A RECOMMENDATION. DECISIONS REGARDING WHETHER AND HOW TO UTILIZE ANY OF THESE PRACTICES SHOULD BE MADE BY HEALTH CARE PROVIDERS, AT THEIR OWN RISK, WITH CONSIDERATION OF INDIVIDUAL CIRCUMSTANCES. AS INFORMATION IS CHANGING RAPIDLY, VIZIENT ENCOURAGES YOU TO ALWAYS REFER TO THE CDC, YOUR STATE'S DEPARTMENT OF HEALTH, AND YOUR LOCAL PUBLIC HEALTH AUTHORITY FOR GUIDANCE. VIZIENT DOES NOT PROVIDE LEGAL, REGULATORY, OR MEDICAL ADVICE AND DISCLAIMS LIABILITY OR RESPONSIBILITY FOR THE ACCURACY, COMPLETENESS, AND/OR CLINICAL EFFICACY AND SAFETY FOR THE PRODUCTS OR PROCESSES CONTAINED HEREIN. MEMBERS SHOULD SEEK THEIR LEGAL COUNSEL'S ADVICE ON LOCAL, STATE, AND FEDERAL LEGAL/REGULATORY MATTERS.

Product overview

Pneumatic tourniquet cuffs are inflatable, single- or dual-bladder cuffs, designed to apply pressure on arterial blood flow in a limb to create a bloodless surgical field. Pneumatic tourniquets provide greater safety to the limb because they can occlude arterial blood flow at lower pressures, thus increasing the probability of saving the limb and recovering full limb.

Tourniquet cuff systems consist of:

- pressure regulated control unit
- hose assembly
- inflatable tourniquet cuff

Best Practices for Safe Use

- [Association of Surgical Technologists \(AST\) Guidelines for Best Practices for Safe Use of Pneumatic Tourniquets](#)
The Certified Surgical Technologists (CST) should know if the patient's surgery involves the use of a pneumatic tourniquet.
 - The surgeon and/or anesthesia provider make the decision if a pneumatic tourniquet is or is not required based upon the results of the patient's history and physical contraindications such as existing disease condition (s) and type and length of surgical procedure.
 - CSTs should not assume that a tourniquet will be used for surgical procedures on the extremities. Studies have provided evidence that the use of a tourniquet in all instances of extremity surgery may be unnecessary, and the CST should confirm with the surgeon and/or anesthesia provider if a tourniquet will be used.
- [Pneumatic Tourniquets \(Association for perioperative Registered Nurses \[AORN\]\)](#)
- AORN Guideline Quick View: Pneumatic Tourniquets (AORN, Inc. 2020)
 - During the preoperative briefing, discuss plans for use of a pneumatic tourniquet
 - Preoperatively assess the patient's
 - Circulation in the extremity (i.e., pulses, temperature, capillary refill)
 - Skin condition under the distal to the planned cuff site
 - Sensory and motor responses of the operative extremity compared to the contralateral extremity
 - Assess the patient preoperatively for factors that may increase his or her risk of developing complications from pneumatic tourniquet use
 - Use pneumatic tourniquets with caution for patients who:
 - Are older adults
 - Have low preoperative hemoglobin
 - Have an increased risk for venous thromboembolism (VTE)
 - Assess the patient for potential systemic contraindications to tourniquet use, including:
 - Arteriovenous grafts or fistulas
 - Diabetic neuropathy
 - Peripheral vascular disease
 - Previous revascularization
 - Sickle cell anemia
 - Severe infection
 - VTE or a history of VTE
 - Notify the prescriber (surgeon or anesthesia professional) of any contraindications to tourniquet use, including:
 - malignancy
 - open fracture
 - severe crushing injuries

Additional resources

[Supply assurance webpage](#); [Vizient Newsroom](#)



Want to receive weekly Supply Assurance updates?

Update your preferences through our [Subscription Manager](#) by selecting Supply Assurance Weekly Digest.

Questions? Contact disasterresponse@vizientinc.com

Competitive cross reference

Stryker non-sterile to competitor sterile cross reference

Stryker PN	Description	Supplier	Part number	Description2
5921-030-135NS	Disposable 30 x 4, single bladder, single port, quick connect	ASP	ASP-564010	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 30X4" 1 PORT (10) WHITE, WITH SLEEVE
	Disposable 30 x 4, single bladder, single port, quick connect	Zimmer	60707510500	30 SPSB STER DISP CUFF PLC
	Disposable 30 x 4, single bladder, single port, quick connect	Zimmer	60707515500	30 SPSB STER DISP CUFF PLC W/SLEEVE
	Disposable 30 x 4, single bladder, single port, quick connect	Delfi	N/A	See sub for 18" or 34"
	Disposable 30 x 4, single bladder, single port, quick connect	Delfi	N/A	See sub for 18" or 34"
	Disposable 30 x 4, single bladder, single port, quick connect	SourceMark	REF 20-34-727SLZ-1	Dispo Cuff, for single use, box of 10
	Disposable 30 x 4, single bladder, single port, quick connect	SourceMark	REF 20-34-727SLZ-1	Dispo Cuff, for single use, box of 10
	Disposable 30 x 4, single bladder, single port, quick connect	SourceMark	REF 20-54-527	Fabric Cuff, washable, box of 1
	Disposable 30 x 4, single bladder, single port, quick connect	SourceMark	REF 20-54-527	Fabric Cuff, washable, box of 1
	Disposable 30 x 4, single bladder, single port, quick connect	SourceMark	REF 20-64-527	Silicone Velcro Cuff, Autoclavable, box of 1
	Disposable 30 x 4, single bladder, single port, quick connect	SourceMark	REF 20-64-527	Silicone Velcro Cuff, Autoclavable, box of 1
	Disposable 30 x 4, single bladder, single port, quick connect	SourceMark	REF 20-74-727	Slim cuff, washable, box of 1
	Disposable 30 x 4, single bladder, single port, quick connect	SourceMark	REF 20-74-727	Slim cuff, washable, box of 1
	Disposable 30 x 4, single bladder, single port, quick connect	SourceMark	REF 20-74-727	Slim cuff, washable, box of 1
5921-030-235NS	Disposable 30 x 4, single bladder, dual port, quick connect	ASP	ASP-564011	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 30X4" 2 PORT (10) WHITE, WITH SLEEVE
	Disposable 30 x 4, single bladder, dual port, quick connect	Zimmer	60707010500	30 DPSB STER DISP CUFF PLC
	Disposable 30 x 4, single bladder, dual port, quick connect	Zimmer	60707015500	30 DPSB STER DISP CUFF PLC W/SLEEVE
	Disposable 30 x 4, single bladder, dual port, quick connect	Delfi	N/A	See sub for 18" or 34"
	Disposable 30 x 4, single bladder, dual port, quick connect	Delfi	N/A	See sub for 18" or 34"
	Disposable 30 x 4, single bladder, dual port, quick connect	SourceMark	REF 20-50-727	Fabric Cuff, washable, box of 1
	Disposable 30 x 4, single bladder, dual port, quick connect	SourceMark	REF 20-50-727	Fabric Cuff, washable, box of 1
5921-034-135NS	Disposable 34 x 4, single bladder, single port, quick connect	ASP	ASP-564012	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 34X4" 1 PORT (10) BLUE, WITH SLEEVE
	Disposable 34 x 4, single bladder, single port, quick connect	Zimmer	60707510600	34 SPSB STER DISP CUFF PLC
	Disposable 34 x 4, single bladder, single port, quick connect	Zimmer	60707515600	34 SPSB STER DISP CUFF PLC W/SLEEVE
	Disposable 34 x 4, single bladder, single port, quick connect	Delfi	9-7500-034	10 per box, Easi-Fit Variable Contour and Cylindrical, Reuseable, 34" x 4.5"
	Disposable 34 x 4, single bladder, single port, quick connect	Delfi	9-7500-034	10 per box, Easi-Fit Variable Contour and Cylindrical, Reuseable, 34" x 4.5"
	Disposable 34 x 4, single bladder, single port, quick connect	SourceMark	REF 20-34-728SLZ-1	Dispo Cuff, for single use, box of 10
	Disposable 34 x 4, single bladder, single port, quick connect	SourceMark	REF 20-34-728SLZ-1	Dispo Cuff, for single use, box of 10
	Disposable 34 x 4, single bladder, single port, quick connect	SourceMark	REF 20-54-528	Fabric Cuff, washable, box of 1
	Disposable 34 x 4, single bladder, single port, quick connect	SourceMark	REF 20-54-528	Fabric Cuff, washable, box of 1
	Disposable 34 x 4, single bladder, single port, quick connect	SourceMark	REF 20-64-528	Silicone Velcro Cuff, Autoclavable, box of 1
	Disposable 34 x 4, single bladder, single port, quick connect	SourceMark	REF 20-64-528	Silicone Velcro Cuff, Autoclavable, box of 1
	Disposable 34 x 4, single bladder, single port, quick connect	SourceMark	REF 20-64-528	Silicone Velcro Cuff, Autoclavable, box of 1

Stryker PN	Description	Supplier	Part number	Description2
	Disposable 34 x 4, single bladder, single port, quick connect	SourceMark	REF 20-74-728	Slim cuff, washable, box of 1
	Disposable 34 x 4, single bladder, single port, quick connect	SourceMark	REF 20-74-728	Slim cuff, washable, box of 1
5921-034-235NS	Disposable 34 x 4, single bladder, dual port, quick connect	ASP	ASP-564013	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 34X4" 2 PORT (10) BLUE, WITH SLEEVE
	Disposable 34 x 4, single bladder, dual port, quick connect	Zimmer	60707010600	34 DPSB STER DISP CUFF PLC
	Disposable 34 x 4, single bladder, dual port, quick connect	Zimmer	60707015600	34 DPSB STER DISP CUFF PLC W/SLEEVE
	Disposable 34 x 4, single bladder, dual port, quick connect	Delfi	9-7600-034	10 per box, Easi-Fit Variable Contour and Cylindrical, Reuseable, 34" x 4.5"
	Disposable 34 x 4, single bladder, dual port, quick connect	Delfi	9-7600-034	10 per box, Easi-Fit Variable Contour and Cylindrical, Reuseable, 34" x 4.5"
	Disposable 34 x 4, single bladder, dual port, quick connect	SourceMark	REF 20-50-728	Fabric Cuff, washable, box of 1
	Disposable 34 x 4, single bladder, dual port, quick connect	SourceMark	REF 20-50-728	Fabric Cuff, washable, box of 1
	Disposable 34 x 4, single bladder, dual port, quick connect	SourceMark	REF 20-50-728	Fabric Cuff, washable, box of 1
5921-218-135NS	Disposable 18x4, single bladder, single port, quick connect	ASP	ASP-564006	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 18X4" 1 PORT (10) RED, WITH SLEEVE
	Disposable 18x4, single bladder, single port, quick connect	Zimmer	60707510300	18 SPSB STER DISP CUFF PLC
	Disposable 18x4, single bladder, single port, quick connect	Zimmer	60707515300	18 SPSB STER DISP CUFF PLC W/SLEEVE
	Disposable 18x4, single bladder, single port, quick connect	Delfi	9-7500-018	10 per box, Easi-Fit Variable Contour and Cylindrical, Reuseable, 18" x 4.5"
	Disposable 18x4, single bladder, single port, quick connect	Delfi	9-7500-018	10 per box, Easi-Fit Variable Contour and Cylindrical, Reuseable, 18" x 4.5"
	Disposable 18x4, single bladder, single port, quick connect	SourceMark	REF 20-34-712SLZ-1 (Arm) / REF 20-34-715SLZ-1 (Arm/leg)	Dispo Cuff, for single use, box of 10
	Disposable 18x4, single bladder, single port, quick connect	SourceMark	REF 20-34-712SLZ-1 (Arm) / REF 20-34-715SLZ-1 (Arm/leg)	Dispo Cuff, for single use, box of 10
	Disposable 18x4, single bladder, single port, quick connect	SourceMark	REF 20-54-712 (Arm) / REF 20-54-512 (Arm/leg)	Fabric Cuff, washable, box of 1
	Disposable 18x4, single bladder, single port, quick connect	SourceMark	REF 20-54-712 (Arm) / REF 20-54-512 (Arm/leg)	Fabric Cuff, washable, box of 1
	Disposable 18x4, single bladder, single port, quick connect	SourceMark	REF 20-64-612 (Arm) / REF 20-64-512 (Arm/leg)	Silicone Velcro Cuff, Autoclavable, box of 1
	Disposable 18x4, single bladder, single port, quick connect	SourceMark	REF 20-64-612 (Arm) / REF 20-64-512 (Arm/leg)	Silicone Velcro Cuff, Autoclavable, box of 1
	Disposable 18x4, single bladder, single port, quick connect	SourceMark	REF 20-74-712 (Arm) / REF 20-74-715 (Arm/leg)	Slim cuff, washable, box of 1
	Disposable 18x4, single bladder, single port, quick connect	SourceMark	REF 20-74-712 (Arm) / REF 20-74-715 (Arm/leg)	Slim cuff, washable, box of 1
5921-218-235NS	Disposable 18x4, single bladder, dual port, quick connect	ASP	ASP-564007	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 18X4" 2 PORT (10) RED, WITH SLEEVE
	Disposable 18x4, single bladder, dual port, quick connect	Zimmer	60707010300	18 DPSB STER DISP CUFF PLC
	Disposable 18x4, single bladder, dual port, quick connect	Zimmer	60707015300	18 DPSB STER DISP CUFF PLC W/SLEEVE
	Disposable 18x4, single bladder, dual port, quick connect	Delfi	9-7600-018	10 per box, Easi-Fit Variable Contour and Cylindrical, Reuseable, 18" x 4.5"
	Disposable 18x4, single bladder, dual port, quick connect	Delfi	9-7600-018	10 per box, Easi-Fit Variable Contour and Cylindrical, Reuseable, 18" x 4.5"
	Disposable 18x4, single bladder, dual port, quick connect	SourceMark	REF 20-30-712SLZ	Dispo Cuff, for single use, box of 10
	Disposable 18x4, single bladder, dual port, quick connect	SourceMark	REF 20-30-712SLZ	Dispo Cuff, for single use, box of 10
	Disposable 18x4, single bladder, dual port, quick connect	SourceMark	REF 20-50-712	Fabric Cuff, washable, box of 1
	Disposable 18x4, single bladder, dual port, quick connect	SourceMark	REF 20-50-712	Fabric Cuff, washable, box of 1
	Disposable 18x4, single bladder, dual port, quick connect	SourceMark	REF 20-60-712	Silicone Velcro Cuff, Autoclavable, box of 1
	Disposable 18x4, single bladder, dual port, quick connect	SourceMark	REF 20-60-712	Silicone Velcro Cuff, Autoclavable, box of 1
	Disposable 18x4, single bladder, dual port, quick connect	SourceMark	REF 20-60-712	Silicone Velcro Cuff, Autoclavable, box of 1

Stryker sterile to competitor sterile cross reference

Stryker PN	Description	Supplier	Part number	Description2
5921-018-145	DISP 18X5.5,2BLA,1PRT QUICK	Zimmer	60708500200	18 IN SPDB STER DISP CUFF LUER LOCK
5921-018-235	DISP 18X3,1BLA,2PRT QUICK	SourceMark	20-36-712SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Long Arm, 18X4 inches, 46X10 cm
5921-018-245	DISP 18X5.5,2BLA,2PRT QUICK	Zimmer	60708000200	18 IN DPDB STER DISP CUFF LUER LOCK
	DISP 18X5.5,2BLA,2PRT QUICK	Zimmer	60708010200	18 IN DPDB STER DISP CUFF PLC
	DISP 18X5.5,2BLA,2PRT QUICK	Zimmer	60708010200	18 IN DPDB STER DISP CUFF PLC
5921-024-135	DISP 24X4,1BLA,1PRT QUICK	ASP	ASP-564008	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 24X4" 1 PORT (10) YELLOW, WITH SLEEVE
	DISP 24X4,1BLA,1PRT QUICK	ASP	ASP-564008	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 24X4" 1 PORT (10) YELLOW, WITH SLEEVE
	DISP 24X4,1BLA,1PRT QUICK	ASP	ASP-564008	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 24X4" 1 PORT (10) YELLOW, WITH SLEEVE
	DISP 24X4,1BLA,1PRT QUICK	ASP	ASP-564008	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 24X4" 1 PORT (10) YELLOW, WITH SLEEVE
	DISP 24X4,1BLA,1PRT QUICK	ASP	ASP-564008	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 24X4" 1 PORT (10) YELLOW, WITH SLEEVE
	DISP 24X4,1BLA,1PRT QUICK	ASP	ASP-564008	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 24X4" 1 PORT (10) YELLOW, WITH SLEEVE
	DISP 24X4,1BLA,1PRT QUICK	SourceMark	20-34-722SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Leg, Conical, 24X4 3/4 inches, 61X12 cm
	DISP 24X4,1BLA,1PRT QUICK	SourceMark	20-34-722SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Leg, Conical, 24X4 3/4 inches, 61X12 cm
	DISP 24X4,1BLA,1PRT QUICK	Zimmer	60707500400	24 IN SPSB STER DISP CUFF LUER LOCK
	DISP 24X4,1BLA,1PRT QUICK	Zimmer	60707510400	24 IN SPSB STER DISP CUFF PLC
	DISP 24X4,1BLA,1PRT QUICK	Zimmer	60707515400	24 IN SPSB STER DISP CUFF PLC W/SLEEVE
	DISP 24X4,1BLA,1PRT QUICK	Zimmer	60707510400	24 IN SPSB STER DISP CUFF PLC
	DISP 24X4,1BLA,1PRT QUICK	Zimmer	60707515400	24 IN SPSB STER DISP CUFF PLC W/SLEEVE
	DISP 24X4,1BLA,1PRT QUICK	Zimmer	60707500400	24 IN SPSB STER DISP CUFF LUER LOCK
5921-024-145	DISP 24X5.5,2BLA,1PRT QUICK	SourceMark	20-30-722SLZ-1	Disposable Cuff, Double Cuffs, Single Port, Leg, 24X6 inches, 61X15 cm
	DISP 24X5.5,2BLA,1PRT QUICK	Zimmer	60708500300	24 IN SPDB STER DISP CUFF LUER LOCK
	DISP 24X5.5,2BLA,1PRT QUICK	Zimmer	60708510300	24 IN SPDB STER DISP CUFF PLC
5921-024-235	DIS 24"X4",1BLA,2PRT QUICK	ASP	ASP-564009	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 24X4" 2 PORT (10) YELLOW, WITH SLEEVE
	DIS 24"X4",1BLA,2PRT QUICK	ASP	ASP-564009	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 24X4" 2 PORT (10) YELLOW, WITH SLEEVE
	DIS 24"X4",1BLA,2PRT QUICK	ASP	ASP-564009	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 24X4" 2 PORT (10) YELLOW, WITH SLEEVE
	DIS 24"X4",1BLA,2PRT QUICK	ASP	ASP-564009	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 24X4" 2 PORT (10) YELLOW, WITH SLEEVE
	DIS 24"X4",1BLA,2PRT QUICK	SourceMark	20-36-722SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Leg, Conical, 24X4 3/4 inches, 61X12 cm
	DIS 24"X4",1BLA,2PRT QUICK	SourceMark	20-36-722SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Leg, Conical, 24X4 3/4 inches, 61X12 cm
	DIS 24"X4",1BLA,2PRT QUICK	Zimmer	60707000400	24 IN DPSB STER DISP CUFF LUER LOCK
	DIS 24"X4",1BLA,2PRT QUICK	Zimmer	60707010400	24 IN DPSB STER DISP CUFF PLC
	DIS 24"X4",1BLA,2PRT QUICK	Zimmer	60707015400	24 IN DPSB STER DISP CUFF PLC W/SLEEVE
	DIS 24"X4",1BLA,2PRT QUICK	Zimmer	60707010400	24 IN DPSB STER DISP CUFF PLC
	DIS 24"X4",1BLA,2PRT QUICK	Zimmer	60707015400	24 IN DPSB STER DISP CUFF PLC W/SLEEVE
5921-024-245	DISP 24X5.5,2BLA,2PRT QUICK	ASP	ASP-564020	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 26X6" DUAL BLADER 2 PORT (10) ORANGE, WITH SLEEVE
	DISP 24X5.5,2BLA,2PRT QUICK	ASP	ASP-564020	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 26X6" DUAL BLADER 2 PORT (10) ORANGE, WITH SLEEVE
	DISP 24X5.5,2BLA,2PRT QUICK	Zimmer	60708000300	24 IN DPDB STER DISP CUFF LUER LOCK
	DISP 24X5.5,2BLA,2PRT QUICK	Zimmer	60708010300	24 IN DPDB STER DISP CUFF PLC

Stryker PN	Description	Supplier	Part number	Description2
	DISP 24X5.5,2BLA,2PRT QUICK	Zimmer	60708010300	24 IN DPDB STER DISP CUFF PLC
5921-030-135	DISP 30X4,1BLA,1PRT QUICK	ASP	ASP-564010	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 30X4" 1 PORT (10) WHITE, WITH SLEEVE
	DISP 30X4,1BLA,1PRT QUICK	ASP	ASP-564010	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 30X4" 1 PORT (10) WHITE, WITH SLEEVE
	DISP 30X4,1BLA,1PRT QUICK	ASP	ASP-564010	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 30X4" 1 PORT (10) WHITE, WITH SLEEVE
	DISP 30X4,1BLA,1PRT QUICK	ASP	ASP-564010	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 30X4" 1 PORT (10) WHITE, WITH SLEEVE
	DISP 30X4,1BLA,1PRT QUICK	ASP	ASP-564010	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 30X4" 1 PORT (10) WHITE, WITH SLEEVE
	DISP 30X4,1BLA,1PRT QUICK	SourceMark	20-34-727SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Long Leg, Conical, 30X4 3/4 inches, 76X12 cm
	DISP 30X4,1BLA,1PRT QUICK	SourceMark	20-34-727SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Long Leg, Conical, 30X4 3/4 inches, 76X12 cm
	DISP 30X4,1BLA,1PRT QUICK	Zimmer	60707500500	30 IN SPSB STER DISP CUFF LUER LOCK
	DISP 30X4,1BLA,1PRT QUICK	Zimmer	60707510500	30 IN SPSB STER DISP CUFF PLC
	DISP 30X4,1BLA,1PRT QUICK	Zimmer	60707515500	30 IN SPSB STER DISP CUFF PLC W/SLEEVE
	DISP 30X4,1BLA,1PRT QUICK	Zimmer	60707510500	30 IN SPSB STER DISP CUFF PLC
	DISP 30X4,1BLA,1PRT QUICK	Zimmer	60707515500	30 IN SPSB STER DISP CUFF PLC W/SLEEVE
	DISP 30X4,1BLA,1PRT QUICK	Zimmer	60707500500	30 IN SPSB STER DISP CUFF LUER LOCK
5921-030-235	DISP 30X4,1BLA,2PRT QUICK	ASP	ASP-564021	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 32X6" DUAL BLADER 1 PORT (10) BROWN, WITH SLEEVE
	DISP 30X4,1BLA,2PRT QUICK	ASP	ASP-564011	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 30X4" 2 PORT (10) WHITE, WITH SLEEVE
	DISP 30X4,1BLA,2PRT QUICK	ASP	ASP-564011	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 30X4" 2 PORT (10) WHITE, WITH SLEEVE
	DISP 30X4,1BLA,2PRT QUICK	ASP	ASP-564011	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 30X4" 2 PORT (10) WHITE, WITH SLEEVE
	DISP 30X4,1BLA,2PRT QUICK	ASP	ASP-564011	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 30X4" 2 PORT (10) WHITE, WITH SLEEVE
	DISP 30X4,1BLA,2PRT QUICK	SourceMark	20-36-727SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Long Leg, Conical, 30X4 3/4 inches, 76X12 cm
	DISP 30X4,1BLA,2PRT QUICK	SourceMark	20-36-727SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Long Leg, Conical, 30X4 3/4 inches, 76X12 cm
	DISP 30X4,1BLA,2PRT QUICK	Zimmer	60707000500	30 IN DPSB STER DISP CUFF LUER LOCK
	DISP 30X4,1BLA,2PRT QUICK	Zimmer	60707010500	30 IN DPSB STER DISP CUFF PLC
	DISP 30X4,1BLA,2PRT QUICK	Zimmer	60707015500	30 IN DPSB STER DISP CUFF PLC W/SLEEVE
	DISP 30X4,1BLA,2PRT QUICK	Zimmer	60707010500	30 IN DPSB STER DISP CUFF PLC
	DISP 30X4,1BLA,2PRT QUICK	Zimmer	60707015500	30 IN DPSB STER DISP CUFF PLC W/SLEEVE
	DISP 30X4,1BLA,2PRT QUICK	Zimmer	60707015500	30 IN DPSB STER DISP CUFF PLC W/SLEEVE
5921-034-135	DISP 34X4,1BLA,1PRT QUICK	ASP	ASP-564012	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 34X4" 1 PORT (10) BLUE, WITH SLEEVE
	DISP 34X4,1BLA,1PRT QUICK	ASP	ASP-564012	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 34X4" 1 PORT (10) BLUE, WITH SLEEVE
	DISP 34X4,1BLA,1PRT QUICK	ASP	ASP-564012	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 34X4" 1 PORT (10) BLUE, WITH SLEEVE
	DISP 34X4,1BLA,1PRT QUICK	ASP	ASP-564012	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 34X4" 1 PORT (10) BLUE, WITH SLEEVE
	DISP 34X4,1BLA,1PRT QUICK	ASP	ASP-564012	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 34X4" 1 PORT (10) BLUE, WITH SLEEVE
	DISP 34X4,1BLA,1PRT QUICK	SourceMark	20-34-728SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Extra Long Leg, Conical, 34X4 3/4 inches, 86X12 cm
	DISP 34X4,1BLA,1PRT QUICK	SourceMark	20-34-728SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Extra Long Leg, Conical, 34X4 3/4 inches, 86X12 cm
	DISP 34X4,1BLA,1PRT QUICK	Zimmer	60707500600	34 IN SPSB STER DISP CUFF LUER LOCK
	DISP 34X4,1BLA,1PRT QUICK	Zimmer	60707510600	34 IN SPSB STER DISP CUFF PLC
	DISP 34X4,1BLA,1PRT QUICK	Zimmer	60707515600	34 IN SPSB STER DISP CUFF PLC W/SLEEVE
	DISP 34X4,1BLA,1PRT QUICK	Zimmer	60707510600	34 IN SPSB STER DISP CUFF PLC
	DISP 34X4,1BLA,1PRT QUICK	Zimmer	60707515600	34 IN SPSB STER DISP CUFF PLC W/SLEEVE
	DISP 34X4,1BLA,1PRT QUICK	Zimmer	60707515600	34 IN SPSB STER DISP CUFF PLC W/SLEEVE

Stryker PN	Description	Supplier	Part number	Description2
5921-034-235	DISP 34X4,1BLA,2PRT QUICK	ASP	ASP-564013	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 34X4" 2 PORT (10) BLUE, WITH SLEEVE
	DISP 34X4,1BLA,2PRT QUICK	ASP	ASP-564013	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 34X4" 2 PORT (10) BLUE, WITH SLEEVE
	DISP 34X4,1BLA,2PRT QUICK	ASP	ASP-564013	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 34X4" 2 PORT (10) BLUE, WITH SLEEVE
	DISP 34X4,1BLA,2PRT QUICK	ASP	ASP-564013	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 34X4" 2 PORT (10) BLUE, WITH SLEEVE
	DISP 34X4,1BLA,2PRT QUICK	ASP	ASP-564013	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 34X4" 2 PORT (10) BLUE, WITH SLEEVE
	DISP 34X4,1BLA,2PRT QUICK	SourceMark	20-36-728SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Conical, Extra Long Leg, 34X4 3/4 inches, 86X12 cm
	DISP 34X4,1BLA,2PRT QUICK	SourceMark	20-36-728SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Conical, Extra Long Leg, 34X4 3/4 inches, 86X12 cm
	DISP 34X4,1BLA,2PRT QUICK	Zimmer	60707000600	34 IN DPSB STER DISP CUFF LUER LOCK
	DISP 34X4,1BLA,2PRT QUICK	Zimmer	60707010600	34 IN DPSB STER DISP CUFF PLC
	DISP 34X4,1BLA,2PRT QUICK	Zimmer	60707015600	34 IN DPSB STER DISP CUFF PLC W/SLEEVE
	DISP 34X4,1BLA,2PRT QUICK	Zimmer	60707010600	34 IN DPSB STER DISP CUFF PLC
	DISP 34X4,1BLA,2PRT QUICK	Zimmer	60707015600	34 IN DPSB STER DISP CUFF PLC W/SLEEVE
5921-044-135	DISP 44X4,1BLA,1PRT QUICK	ASP	ASP-564014	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 44X4" 1 PORT (10) GREEN, WITH SLEEVE
	DISP 44X4,1BLA,1PRT QUICK	ASP	ASP-564014	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 44X4" 1 PORT (10) GREEN, WITH SLEEVE
	DISP 44X4,1BLA,1PRT QUICK	ASP	ASP-564014	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 44X4" 1 PORT (10) GREEN, WITH SLEEVE
	DISP 44X4,1BLA,1PRT QUICK	SourceMark	20-34-729SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Super Long Leg, Conical, 42X4 3/4 inches, 107X12 cm
	DISP 44X4,1BLA,1PRT QUICK	Zimmer	60707500700	42 IN SPSB STER DISP CUFF LUER LOCK
	DISP 44X4,1BLA,1PRT QUICK	Zimmer	60707510700	42 IN SPSB STER DISP CUFF PLC
	DISP 44X4,1BLA,1PRT QUICK	Zimmer	60707515700	42 IN SPSB STER DISP CUFF PLC W/SLEEVE
	DISP 44X4,1BLA,1PRT QUICK	Zimmer	60707515700	42 IN SPSB STER DISP CUFF PLC W/SLEEVE
5921-044-235	DISP 44X4,1BLA,2PRT QUICK	ASP	ASP-564015	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 44X4" 2 PORT (10) GREEN, WITH SLEEVE
	DISP 44X4,1BLA,2PRT QUICK	ASP	ASP-564015	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 44X4" 2 PORT (10) GREEN, WITH SLEEVE
	DISP 44X4,1BLA,2PRT QUICK	ASP	ASP-564015	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 44X4" 2 PORT (10) GREEN, WITH SLEEVE
	DISP 44X4,1BLA,2PRT QUICK	SourceMark	20-36-729SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Super Long Leg, Conical, 42X4 3/4 inches, 107X12 cm
	DISP 44X4,1BLA,2PRT QUICK	SourceMark	20-36-729SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Super Long Leg, Conical, 42X4 3/4 inches, 107X12 cm
	DISP 44X4,1BLA,2PRT QUICK	Zimmer	60707000700	42 IN DPSB STER DISP CUFF LUER LOCK
	DISP 44X4,1BLA,2PRT QUICK	Zimmer	60707010700	42 IN DPSB STER DISP CUFF PLC
	DISP 44X4,1BLA,2PRT QUICK	Zimmer	60707015700	42 IN DPSB STER DISP CUFF PLC W/SLEEVE
	DISP 44X4,1BLA,2PRT QUICK	Zimmer	60707010700	42 IN DPSB STER DISP CUFF PLC
	DISP 44X4,1BLA,2PRT QUICK	Zimmer	60707015700	42 IN DPSB STER DISP CUFF PLC W/SLEEVE
	DISP 44X4,1BLA,2PRT QUICK	Zimmer	60707015700	42 IN DPSB STER DISP CUFF PLC W/SLEEVE
5921-115-135	DISP 15X2.5,1BLA,1PRT QUICK	ASP	ASP-564004	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 15X4" 1 PORT (10) GRAY, WITH SLEEVE
	DISP 15X2.5,1BLA,1PRT QUICK	ASP	ASP-564004	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 15X4" 1 PORT (10) GRAY, WITH SLEEVE
	DISP 15X2.5,1BLA,1PRT QUICK	ASP	ASP-564004	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 15X4" 1 PORT (10) GRAY, WITH SLEEVE
	DISP 15X2.5,1BLA,1PRT QUICK	SourceMark	20-34-711SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Arm, 14X3 1/8 inches, 35X8 cm
5921-115-235	DISP 15X2.5,1BLA,2PRT QUICK	ASP	ASP-564005	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 15X4" 2 PORT (10) GRAY, WITH SLEEVE
	DISP 15X2.5,1BLA,2PRT QUICK	ASP	ASP-564005	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 15X4" 2 PORT (10) GRAY, WITH SLEEVE
	DISP 15X2.5,1BLA,2PRT QUICK	ASP	ASP-564005	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 15X4" 2 PORT (10) GRAY, WITH SLEEVE

Stryker PN	Description	Supplier	Part number	Description2
	DISP 15X2.5,1BLA,2PRT QUICK	SourceMark	20-36-711SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Arm, 14X3 1/8 inches, 35X8 cm
	DISP 15X2.5,1BLA,2PRT QUICK	SourceMark	20-36-711SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Arm, 14X3 1/8 inches, 35X8 cm
5921-209-135	DISP 9X2.75,1BLA,1PRT QUICK	ASP	ASP-564000	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 8 X2.5" 1 PORT (10) LIGHT BLUE, WITH SLEEVE
	DISP 9X2.75,1BLA,1PRT QUICK	ASP	ASP-564000	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 8 X2.5" 1 PORT (10) LIGHT BLUE, WITH SLEEVE
	DISP 9X2.75,1BLA,1PRT QUICK	SourceMark	20-34-700SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Infant, 8X2 inches, 20X5 cm
	DISP 9X2.75,1BLA,1PRT QUICK	SourceMark	20-34-700SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Infant, 8X2 inches, 20X5 cm
	DISP 9X2.75,1BLA,1PRT QUICK	SourceMark	20-34-700SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Infant, 8X2 inches, 20X5 cm
	DISP 9X2.75,1BLA,1PRT QUICK	Zimmer	60707500100	8 IN SPSB STER DISP CUFF LUER LOCK
	DISP 9X2.75,1BLA,1PRT QUICK	Zimmer	60707510100	8 IN SPSB STER DISP CUFF PLC
	DISP 9X2.75,1BLA,1PRT QUICK	Zimmer	60707510100	8 IN SPSB STER DISP CUFF PLC
5921-209-235	DISP 9X2.75,1BLA,2PRT QUICK	SourceMark	20-36-700SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Infant, 8X2 inches, 20X5 cm
	DISP 9X2.75,1BLA,2PRT QUICK	SourceMark	20-36-700SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Infant, 8X2 inches, 20X5 cm
	DISP 9X2.75,1BLA,2PRT QUICK	SourceMark	20-36-700SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Infant, 8X2 inches, 20X5 cm
	DISP 9X2.75,1BLA,2PRT QUICK	Zimmer	60707000100	8 IN DPSB STER DISP CUFF LUER LOCK
	DISP 9X2.75,1BLA,2PRT QUICK	Zimmer	60707010100	8 IN DPSB STER DISP CUFF PLC
	DISP 9X2.75,1BLA,2PRT QUICK	Zimmer	60707015100	8 DPSB STER DISP CUFF PLC W/SLEEVE
	DISP 9X2.75,1BLA,2PRT QUICK	Zimmer	60707010100	8 IN DPSB STER DISP CUFF PLC
	DISP 9X2.75,1BLA,2PRT QUICK	Zimmer	60707010100	8 IN DPSB STER DISP CUFF PLC
5921-212-135	DISP 12X3.5,1BLA,1PRT QUICK	ASP	ASP-564002	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 12X3.5" 1 PORT (10) PLUM, WITH SLEEVE
	DISP 12X3.5,1BLA,1PRT QUICK	ASP	ASP-564002	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 12X3.5" 1 PORT (10) PLUM, WITH SLEEVE
	DISP 12X3.5,1BLA,1PRT QUICK	ASP	ASP-564002	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 12X3.5" 1 PORT (10) PLUM, WITH SLEEVE
	DISP 12X3.5,1BLA,1PRT QUICK	ASP	ASP-564002	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 12X3.5" 1 PORT (10) PLUM, WITH SLEEVE
	DISP 12X3.5,1BLA,1PRT QUICK	ASP	ASP-564002	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 12X3.5" 1 PORT (10) PLUM, WITH SLEEVE
	DISP 12X3.5,1BLA,1PRT QUICK	SourceMark	20-34-710SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Child, 12X2 3/8 inches, 30X6 cm
	DISP 12X3.5,1BLA,1PRT QUICK	SourceMark	20-34-710SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Child, 12X2 3/8 inches, 30X6 cm
	DISP 12X3.5,1BLA,1PRT QUICK	SourceMark	20-34-710SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Child, 12X2 3/8 inches, 30X6 cm
	DISP 12X3.5,1BLA,1PRT QUICK	SourceMark	20-34-710SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Child, 12X2 3/8 inches, 30X6 cm
	DISP 12X3.5,1BLA,1PRT QUICK	Zimmer	60707500200	12 IN SPSB STER DISP CUFF LUER LOCK
	DISP 12X3.5,1BLA,1PRT QUICK	Zimmer	60707510200	12 IN SPSB STER DISP CUFF PLC
	DISP 12X3.5,1BLA,1PRT QUICK	Zimmer	60707515200	12 IN SPSB STER DISP CUFF PLC W/SLEEVE
	DISP 12X3.5,1BLA,1PRT QUICK	Zimmer	60707510200	12 IN SPSB STER DISP CUFF PLC
	DISP 12X3.5,1BLA,1PRT QUICK	Zimmer	60707515200	12 IN SPSB STER DISP CUFF PLC W/SLEEVE
	DISP 12X3.5,1BLA,1PRT QUICK	Zimmer	60707510200	12 IN SPSB STER DISP CUFF PLC
5921-212-235	DISP 12X3.5,1BLA,2PRT QUICK	ASP	ASP-564003	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 12X3.5" 2 PORT (10) PLUM, WITH SLEEVE
	DISP 12X3.5,1BLA,2PRT QUICK	ASP	ASP-564003	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 12X3.5" 2 PORT (10) PLUM, WITH SLEEVE
	DISP 12X3.5,1BLA,2PRT QUICK	ASP	ASP-564003	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 12X3.5" 2 PORT (10) PLUM, WITH SLEEVE
	DISP 12X3.5,1BLA,2PRT QUICK	SourceMark	20-36-710SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Child, 12X2 3/8 inches, 30X6 cm
	DISP 12X3.5,1BLA,2PRT QUICK	SourceMark	20-36-710SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Child, 12X2 3/8 inches, 30X6 cm

Stryker PN	Description	Supplier	Part number	Description2
	DISP 12X3.5,1BLA,2PRT QUICK	SourceMark	20-36-710SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Child, 12X2 3/8 inches, 30X6 cm
	DISP 12X3.5,1BLA,2PRT QUICK	Zimmer	60707000200	12 IN DPSB STER DISP CUFF LUER LOCK
	DISP 12X3.5,1BLA,2PRT QUICK	Zimmer	60707010200	12 IN DPSB STER DISP CUFF PLC
	DISP 12X3.5,1BLA,2PRT QUICK	Zimmer	60707015200	12 IN DPSB STER DISP CUFF PLC W/SLEEVE
	DISP 12X3.5,1BLA,2PRT QUICK	Zimmer	60707010200	12 IN DPSB STER DISP CUFF PLC
5921-218-135	DISP 18X4,1 BLA,1PRT QUICK	ASP	ASP-564006	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 18X4" 1 PORT (10) RED, WITH SLEEVE
	DISP 18X4,1 BLA,1PRT QUICK	ASP	ASP-564006	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 18X4" 1 PORT (10) RED, WITH SLEEVE
	DISP 18X4,1 BLA,1PRT QUICK	ASP	ASP-564006	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 18X4" 1 PORT (10) RED, WITH SLEEVE
	DISP 18X4,1 BLA,1PRT QUICK	ASP	ASP-564004	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 15X4" 1 PORT (10) GRAY, WITH SLEEVE
	DISP 18X4,1 BLA,1PRT QUICK	ASP	ASP-564006	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 18X4" 1 PORT (10) RED, WITH SLEEVE
	DISP 18X4,1 BLA,1PRT QUICK	ASP	ASP-564006	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 18X4" 1 PORT (10) RED, WITH SLEEVE
	DISP 18X4,1 BLA,1PRT QUICK	SourceMark	20-30-712SLZ-1	Disposable Cuff, Double Cuffs, Single Port, Long Arm, 18X6 inches, 46X15 cm
	DISP 18X4,1 BLA,1PRT QUICK	SourceMark	20-34-712SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Long Arm, 18X4 inches, 46X10 cm
	DISP 18X4,1 BLA,1PRT QUICK	SourceMark	20-34-711SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Arm, 14X3 1/8 inches, 35X8 cm
	DISP 18X4,1 BLA,1PRT QUICK	SourceMark	20-34-712SLZ-1	Disposable Cuff, Single Cuff, Single Port, Sterile, Long Arm, 18X4 inches, 46X10 cm
	DISP 18X4,1 BLA,1PRT QUICK	Zimmer	60708510200	18 IN SPDB STER DISP CUFF PLC
	DISP 18X4,1 BLA,1PRT QUICK	Zimmer	60707500300	18 IN SPSB STER DISP CUFF LUER LOCK
	DISP 18X4,1 BLA,1PRT QUICK	Zimmer	60707510300	18 IN SPSB STER DISP CUFF PLC
	DISP 18X4,1 BLA,1PRT QUICK	Zimmer	60707515300	18 IN SPSB STER DISP CUFF PLC W/SLEEVE
	DISP 18X4,1 BLA,1PRT QUICK	Zimmer	60707510300	18 IN SPSB STER DISP CUFF PLC
5921-218-235	DISP 18X4,1BLA,2PRT QUICK	ASP	ASP-564007	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 18X4" 2 PORT (10) RED, WITH SLEEVE
	DISP 18X4,1BLA,2PRT QUICK	ASP	ASP-564007	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 18X4" 2 PORT (10) RED, WITH SLEEVE
	DISP 18X4,1BLA,2PRT QUICK	ASP	ASP-564007	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 18X4" 2 PORT (10) RED, WITH SLEEVE
	DISP 18X4,1BLA,2PRT QUICK	ASP	ASP-564007	TOURNIQUET DISPOSABLE TOURNIQUET CUFF 18X4" 2 PORT (10) RED, WITH SLEEVE
	DISP 18X4,1BLA,2PRT QUICK	SourceMark	20-36-712SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Long Arm, 18X4 inches, 46X10 cm
	DISP 18X4,1BLA,2PRT QUICK	SourceMark	20-36-712SLZ-1	Disposable Cuff, Single Cuff, Dual Port, Sterile, Long Arm, 18X4 inches, 46X10 cm
	DISP 18X4,1BLA,2PRT QUICK	Zimmer	60707000300	18 IN DPSB STER DISP CUFF LUER LOCK
	DISP 18X4,1BLA,2PRT QUICK	Zimmer	60707010300	18 IN DPSB STER DISP CUFF PLC
	DISP 18X4,1BLA,2PRT QUICK	Zimmer	60707015300	18 IN DPSB STER DISP CUFF PLC W/SLEEVE
	DISP 18X4,1BLA,2PRT QUICK	Zimmer	60707010300	18 IN DPSB STER DISP CUFF PLC
5922-008-135	REUSE 8X2.5, 1BLA, 1PRT, QUIK	SourceMark	20-64-700	Reusable Cuff, Silicone/Velcro Single Cuff, Single Port, Infant, 8X1 inches, 20X2.5 cm
5922-012-135	REUSE 12X3, 1BLA, 1PRT, QUIK	SourceMark	20-64-710	Reusable Cuff, Silicone/Velcro Single Cuff, Single Port, Child, 12X1 5/8 inches, 30X4 cm
5922-018-135	REUSE 18X3, 1BLA, 1 PRT, QUICK	SourceMark	20-64-612	Reusable Cuff, Silicone/Velcro Single Cuff, Single Port, Long Arm, 18X3 inches, 46X7.5 cm
5922-018-145	REUSE 18X5.5,2BLA,1PRT,QUIK	SourceMark	20-60-712	Reusable Cuff, Silicone/Velcro Double Cuffs, Long Arm, 18X4 3/4 inches, 46X12 cm
	REUSE 18X5.5,2BLA,1PRT,QUIK	SourceMark	20-60-712	Reusable Cuff, Silicone/Velcro Double Cuffs, Long Arm, 18X4 3/4 inches, 46X12 cm
	REUSE 18X5.5,2BLA,1PRT,QUIK	SourceMark	20-60-712	Reusable Cuff, Silicone/Velcro Double Cuffs, Long Arm, 18X4 3/4 inches, 46X12 cm
5922-024-145	REUSE 24X5.5,2BLA,1PRT,QUIK	SourceMark	20-60-722	Reusable Cuff, Silicone/Velcro Double Cuffs, Leg, 24X5 1/4 inches, 61X13.5 cm
	REUSE 24X5.5,2BLA,1PRT,QUIK	SourceMark	20-60-722	Reusable Cuff, Silicone/Velcro Double Cuffs, Leg, 24X5 1/4 inches, 61X13.5 cm